



Behind Every Smile is a Story


- QUALITY - ЯКІСТЬ
- PRECISION - ТОЧНІСТЬ
- ACCURACY - ПРАВИЛЬНІСТЬ
- DURABILITY - ЗНОСОСТІЙКІСТЬ



BUSHIDO
INSTRUMENTS

- Super-sharp diamond-bead cutter
- smooth ultra-polished and rounded edges
- brilliant stainless titanium stainless steel finish
- superior T-12 stainless steel tool inserts
- sliding lock joint for balance & precision alignment

- Super-sharp diamond-bead cutter
- smooth ultra-polished and rounded edges
- brilliant stainless titanium stainless steel finish
- superior T-12 stainless steel tool inserts
- sliding lock joint for balance & precision alignment



STYLE
A



STYLE
B

ORTHODONTIC INSTRUMENTS

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- أدوات تقويم الأسنان



DE-1800
Height Gauge 0.18
Aluminum



DE-1801
Height Gauge 0.22
Aluminum



DE-1802
Height Gauge 0.18
Stainless Steel



DE-1803
Height Gauge 0.22
Stainless Steel



DE-1804
Moveable
Height Gauge



DE-1805
Endo Finger
Ruler Span
Measure Scale



DE-1806
Endo Finger
Ruler Span
Measure Scale



DE-1807
Endo Finger
Ruler Span
Measure Scale



DE-1808
Positive
Postioner



DE-1809
Strip Saw



DE-1810
Bone Bracket Gauge



DE-1811
Endo Finger Ruler
Span Measure Scale



DE-1812
Endo Finger Ruler
Span Measure Scale

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- ادوات تقويم الأسنان



DE-1813
Measuring Loop
angled, length
175 mm,
measuring range
0 - 40 mm



DE-1814
Measuring Loop
straight length
175 mm,
measuring range
0 - 40 mm



DE-1815
Measuring loop, Castroviejo
angled, length 175 mm,
measuring range 0 - 40 mm,
with window



DE-1816
Measuring Loop
straight length 90 mm,
measuring range
0 - 40 mm,
DE-1817
Castroviejo Calipers
20mm
90mm



DE-1818
Measuring Loop
angled, length 90 mm,
measuring range
0 - 40 mm,
DE-1819
Castroviejo Calipers
20mm, 45°
90mm



DE-1820
Castroviejo Calipers
20mm
90mm



DE-1821
Fig. 1
Iwanson
spring caliper
For metal and
porcelain.



DE-1822
Fig. 2
Iwanson
spring caliper
For wax.

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- أدوات تقويم الأسنان



DE-1823
Boley Gauge



DE-1824
Vernier Gauge



DE-1825
Aluminium Caliper
up to 200mm



DE-1826
Boley Gauge



DE-1827
Delux Boley
Gauge



DE-1828
Beerendonk
130mm
Mini Gauge



DE-1829
Compasso
De Willis



DE-1830
Dentaurum
Ruler



DE-1831
Ruler Wedge



DE-1832
Ruler

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- أدوات تقويم الأسنان

Bracket Holder Tweezers



DE-1833
Deluxe Bracket
Holder Tweezers



DE-1834
Deluxe Bracket
Holder Tweezers



DE-1835
Direct Bond Bracket
Holder Tweezers



DE-1836
Bracket Tweezers
Slot Aligner



DE-1837
Direct Bond
Bracket Holder
Tweezers, Small Tip



DE-1838
Exact Bracket
Holder Tweezers



DE-1839
Bracket Placing
Tweezers



DE-1840
Easy Access
Tweezers

BUSHIDO INSTRUMENTS



DE-1841
Bracket Placer
Tweezers
(Extra Fine Tip)



DE-1842
Deluxe Bracket
Holder Tweezers
135mm



DE-1843
Easy Eccess
Tweezers



DE-1844
Slipper
Tweezers



DE-1845
Direct Bonding
Bracket Holders
135mm



DE-1846
Direct Bonding
Bracket Holders
136mm



DE-1847
Direct Bonding
Bracket Holders
135mm



DE-1848
Deluxe Bracket
Holder



DE-1849
Direct Bonding
Bracket Holders
140mm

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- أدوات تقويم الأسنان



DE-1850
Band seater
auto cleavable



DE-1851
Band seater
auto cleavable
Triangle Shape



DE-1852
Band seater
auto cleavable
Square Shape



DE-1853
Band seater
auto cleavable
Round Shape



DE-1854
Mershon Band
Pusher
140mm



DE-1855
Band seating &
removing
instruments



DE-1856
Ortho Light double
ended
ligature director
scaler



DE-1857
Ortho Light Band
Pusher Scaler



DE-1858
Band Seating &
Removing
Silicon Grip



DE-1859
Ligature Instrument
Silicon Grip



DE-1860
Elastic Instrument
Silicon Grip



DE-1861
Distal Bender
Silicon Grip



DE-1862
Distal Bender
Silicon Grip



DE-1863
Distal Bender
Silicon Grip

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- ادوات تقويم الأسنان



DE-1864
Ortho light scaler double ended



DE-1865
Ortho light scaler double ended



DE-1866
Ortho light scaler double ended



DE-1867
Twiel-on elastic module instrument



DE-1868
Twiel-on elastic module instrument



DE-1869
Schure Band Seater



DE-1870
Ortho light scaler



DE-1871
Ortho light scaler



DE-1872
Ortho light scaler heavy



DE-1873
Ortho instrument



DE-1874
Ortho instrument

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- أدوات تقويم الأسنان



DE-1875
Ortho Light double ended
ligature director



DE-1876
Ortho Light double ended ligature scaler director tivate coil spring
are remove modules with one instruments



DE-1877
Ortho band pusher scaler



DE-1878
Ortho ligature director single ended



DE-1879
Ortho Light double ended ligature
director scaler



DE-1880
Ortho Light double ended ligature scaler/director
tivate coil spring are remove modules with one instruments



DE-1881
Ortho Light double ended ligature director



DE-1882
Ortholite Ligature Director Scaler



DE-1883
Ortho light ligature director double ended



DE-1884
Double ended cement spatula



DE-1885
Ortho light scaler single ended

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتقويم الأسنان برومبیتس

BOX JOINT PLIERS

DE-1886
Universal Dolphin Style
Wire Bending & cutting pliers



DE-1892
Bracket Removing Plier curved
Bracket Removing Plier curved for
easy removal of brackets in the
posterior region. Length 130mm



DE-1887 Band Removing Plier
The curved blade of the band
removing plier can strip easily
and gently remove directly to
the tooth. Length 130mm



DE-1893
Young Wire Bending &
Loop Forming Pliers



DE-1888
Tweed Loop
forming Pliers



DE-1894
Wire Bending &
Loop Forming Pliers



DE-1889
How Plier curved
Universally usable tool for
the multiband technology for
wires up to .028".
Length 130mm



DE-1895
Adams, small wire bending
and loop forming Pliers



DE-1890 How Plier
Universally usable tool for the
multiband technology for wires
up to .028".
Length 135mm



DE-1896
Bracket Removing Plier Round Box
Bracket removing plier for comfortable
and quick removal of the brackets.
Length 135mm



DE-1891
Nance Loop
forming Pliers



DE-1897
Weingart Plier medium
Weingart plier that are ideal due to
their curved tip for setting sheets and
is capable of bending the ends.
Length 145mm



- Щипцы для ортодонтии и протезирования
- Пинцес pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمانشة لتقويم الأسنان بروسیتیس

BOX JOINT PLIERS

DE-1898 Adams Plier Round Box
Adams pliers for bending wire arches and bends right angles for wires up to .022", .025 tip to beak. Length 125mm



DE-1904 Adams Plier
Adams pliers for bending wire arches and bends right angles for wires up to .022", .025 tip to beak. Length 125mm



DE-1899 Falcon Beak Contouring Pliers.



DE-1905 Tweed Loop Formin Pliers



DE-1900
wires up to 0.020" / 0.7mm.



DE-1906 Bird Beak Plier Standard
Bird Beak Wire bending plier without notches, which is suitable as a universal utility tool. Length 125mm



DE-1901 Johnson Contouring Pliers



DE-1907 Mini Ligature cutters



DE-1902 Aderer Large wire bending and loop forming pliers



DE-1908 Small Fluch-cut distal and cutting Pliers
with safety hold fitted with un-breakable, corrosion resistant and durable stainless steel tips



DE-1903 Tweed / angle rectangular arch forming & Torquing pliers



DE-1909 Large Fluch-cut distal and cutting Pliers
with safety hold fitted with un-breakable, corrosion resistant and durable stainless steel tips



- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتقويم الأسنان برومبشيس

BOX JOINT PLIERS

DE-1910
Tweed / Angel
1.78mm (.070") wide tips arch forming and torquing pliers



DE-1916
Contouring Plier



DE-1911
Tweed / Angel
1.78mm (.070") wide tips arch forming and torquing pliers



DE-1917
Bands & Crowns Plier
in order to skirt



DE-1912
Lingual / Palatal arch forming pliers



DE-1918
Ligature Forming Pliers



DE-1913
Aderer
Facebow loop forming Pliers



DE-1919
for double back and triple back bands in 0.30" / 0.76mm & 0.36" / 0.91mm lingual arch wires



DE-1914
Waldsavs wire bending and cutting pliers



DE-1920
Bird Beak Wire Bending & cutting pliers



DE-1915 Aderer Plier
For the manufacturing of appliances and activation of wires for wires up to .030" .fine tip and up to .036" mean tip.
Length 130mm



DE-1921
Weingart Plier medium
Weingart plier that are ideal due to their curved tip for setting sheets and is capable of bending the ends.
Length 145mm



- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتقويم الأسنان بروسیتیس

L-KEY JOINT PLIERS



DE-1922
Distal end cutter with safety hold, mini head, long handles.
fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.64 mm (.018" x .025")
- ⊗ 0.38 mm to 0.53 mm (.0155" to .0215")
- ⊗ Max. 0.53 mm x 0.64 mm (.021" x .025")



DE-1923
Flush-cut distal end cutter with safety hold fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.64 mm (.018" to .025")



DE-1924
Soft wire and ligature cutter, Mini fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ⌋ .25 mm to .38 mm (.010" to .015")



DE-1925
Soft wire and ligature cutter, standard fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ⌋ .25 mm to .38 mm (.010" to .015")



DE-1926
Soft wire and ligature cutter, 15° angled fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ⌋ .25 mm to .38 mm (.010" to .015")



DE-1927
Soft wire and ligature cutter, 45° angled fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ⌋ .25 mm to .38 mm (.010" to .015")



DE-1928
Hard wire cutter, straight Fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



DE-1929
Hard wire cutter, 15° angled Fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



DE-1930
Hard wire cutter, 15° angled Non-autoclavable

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")

- Щипцы для ортодонтии и протезирования
- Щипці для ортодонтії та протезування
- Pincers pour l'orthodontie et la prothèse
- Alicates para ortodoncia y prótesis
- Zangen für Kieferorthopädie und Prothetik
- مكاشة لتقويم الأسنان برومینیس

L-KEY JOINT PLIERS



DE-1931
Band removing pliers



DE-1936
Weingart-mini
Utility Pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")



DE-1932
Bracket removing pliers angled



DE-1937
Weingart
Utility Pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")



DE-1933
Bracket removing pliers straight



DE-1938
How Straight
Utility Pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")



DE-1934
Bracket removing pliers with plastic tip



DE-1939
Angled
Utility Pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")



DE-1935
Adhesive removing pliers



DE-1940
Wire bending and loop forming pliers
Angled (Bird beak)

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتقويم الأسنان بروسیتیس

L-KEY JOINT PLIERS



DE-1941
Multi-purpose pliers, bird beak,
De la rosa and cutter in one

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.56 mm (.018" x .022")



DE-1946
Stop pliers V-Bebd - 1mm

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.43 mm x 0.56 mm (.017" x .028")



DE-1942
Cinch Back / tie back
Distal end bending pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



DE-1947
Lingual / palatal arch forming pliers

- 0.76 mm (.030") / 0.91 mm (.036")



DE-1943
Hammer head
Distal end bending pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



DE-1948
Torquing pliers set



DE-1944
Nance
Loop forming pliers

- 0.36 mm to 0.71 mm (.014" to .028")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



DE-1949
Tweed / angel
Rectangular arch forming and
torquing pliers

- 0.30 mm to 0.71 mm (.012" to .028")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



DE-1945
Detailing step pliers / Bayonet
bending pliers - 1mm

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



DE-1950
Surgical ball hook crimping pliers,
angled



DE-1951
Surgical ball hook crimping
pliers, straight

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمانة لتقويم الأسنان برومینیس


L-KEY JOINT PLIERS




DE-1952
Aderer Clamp three jaw wire bending plier for contouring and bending wire up to .30"/.76mm.



DE-1958
Falcon beak contouring plier.




DE-1953
Bird beak plier for bending round, square and rectangle wire up to .30"/.76mm.



DE-1959
For placement and removal of arch wires, pins and other auxiliaries.




DE-1954
Bird beak plier for light wire bending round, square and rectangular up to .016"/0.41mm.



DE-1960
For placement and removal of arch wires, pins and other auxiliaries in angle design



DE-1955
Double Rounded jaw bending pliers, hooks and wire up to .020"/0.51mm.




DE-1961
contouring Plier.



DE-1962
Bending plier. In order to skirt molares bands. With grooved. Recommended wire. 016", 018" and 022" 0.41mm, 0.46mm and 0.56mm



DE-1956
Adams pliers for hooks and bending wire up to .030"/0.76mm.



DE-1963
Loop forming in four length 3,4,5 & 6 mm wire capacity 0.22", 56mm.



DE-1957
Johnson pliers for crown band contouring



DE-1964
Arch bending plier design for convenient handling of square or rectangular wire, .020" / 0.51mm.



- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتقويم الأسنان بروسیتیس


LAP & L-KEY JOINT PLIERS



DE-1965
In order to skirt bands and crowns.



DE-1966
Wires up to .020" / 0.7mm.



DE-1967
Round and concave plier for loop forming and banding wire up to .20" / .051mm.



DE-1968
Bending for round square or rectangular wire up to 0.30"/0.76mm



DE-1969
for double back and triple back bands in 0.30" / .76mm & 0.36" / .91mm lingual arch wires.



DE-1970
Detailing step plier double-sided offset beak permits bayonet bends. Creates right hand or left hand step-up or stepdown at 1/8, 1/4, or 1mm increments.



DE-1971
"V" bend stop special design makes 1mm "v" bends to shorten archwires or provide a positive stop. Excellent for placing stops in titanium wires. For round or rectangular wires up to .025"



DE-1972
Angle how plier efficient utility plier with serrated beaks for a firm grip. 45 angled tips to access hard to reach areas.



DE-1973
Weingardt plier serrated tips hold wires firmly at a convenient working angle. Tapered beaks fit easily between brackets and hard to reach areas. Our most popular utility plier.



DE-1974
Pin & Ligature Cutters TC.



DE-1975
Mini Pin & Ligature Cutter TC.



DE-1976
Standard Ligature Cutter TC.



DE-1977
V Bend Stop Plier TC.



DE-1978
Angled Direct Bond Bracket Remover TC.



DE-1979
Cutter Str.



Distal Cutter 14cm
DE-1980 2.5 mm
DE-1981 3.0 mm
DE-1982 3.0 mm



DE-1983
Arrow Clasp Forming Pliers 12.0cm



DE-1984
Schwarz 12.5cm



DE-1985
Distal end Cutter, Safety Hold Plier
(Ideal for Cutting Wire) TC.



DE-1986
Hook Crimping Plier Angled
Straight TC.



DE-1987
3 Prong Clasp Adjustable
Plier TC.



DE-1988
Utility Plier Weingart TC.



DE-1989
Direct Bracket Remover
TC.



DE-1990
Straight Direct Bond Bracket
Remover TC.



DE-1991
Pin Holding Plier TC.



DE-1992
Hard wire cutters



DE-1993
Elastic Separating Plier
Round Handle



DE-1994
Elastic
Separating Plier



DE-1995 Str.
DE-1996 Cvd.
Smaha
Elastomeric
ligature
placing forceps
140 mm



DE-1997
Elastomeric
ligature placing
forceps with
groove
Falcon 140 mm



DE-1998
Elastomeric
ligature placing
forceps with
groove
Falcon 140 mm



- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمامة لتقويم الأسنان بروسثیتیس

CUTTER, UTILITY PLIERS & WIRE FORMING PLIERS



DE-1999
Weingart Pliers



DE-2004
Adhesive
Removing Pliers



DE-2000
Slim Weingart
Pliers



DE-2005
Posterior Band Removing
Pliers - Long



DE-2001
How Pliers
Straight



DE-2006
Posterior Band Removing
Pliers - Short




DE-2002
How Pliers Offset



DE-2007
Bird Beak Pliers



DE-2003
Band seating Pliers



Replacement tip and
hex key included.

DE-2008
Omega Loop
Forming Pliers

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонтії та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتقويم الأسنان برومینیس

CUTTERS



DE-2009
Hard
Wire Cutter 15



DE-2013
Mini Pin &
Ligature Cutter



DE-2010
Distal End Cutter
Flush Cut



DE-2014
Hard Wire Cutter
Straight



DE-2011
Universal Cut & Hold
Distal End Cutter



DE-2015
Micro-Mini Pin &
Ligature Cutter



DE-2012
Pin & Ligature
Cutter



DE-2016
Cutter

PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمانشة لتقويم الأسنان بروسیتیس



DE-2017
Three Jaw Pliers



DE-2021
Jarabak Pliers



DE-2018
Light Wore Bird Beak
Pliers with Groove



DE-2022
Hollow Chop
Contouring Pliers



DE-2019
Tweed loop
Forming Pliers



DE-2023
Lingual Arch
Forming Pliers



DE-2020
Arch Bending Pliers



DE-2024
Rectangular Arch
Bending Pliers

ORTHODONTICS PLIERS STANDS

- Ортодонтия Щипцы Подставки
- Supports de pinces d'orthodontie
- Kieferorthopädie Zangenständer
- Ортодонтия Щипці Підставки
- KFO Zangen Ständer
- تقويم الاسنان كمشاة تقف



DE-2025
Plier Stand
Stainless Steel.



DE-2026
Plier Stand
Stainless Steel.



DE-2027
Pliers Stand
Stainless Rod
(Aluminum Anodized)
Stainless Steel



DE-2028
Plier Stand
Stainless Steel.



DE-2029
Plier Stand
Stainless Steel.



DE-2030
Plier Stand
Plastic